 NeoMedLight

# **BiliCocoon™** Fixation System

**Table of contents**

- 1. Product description and application .....5
  - A. Description .....5
  - B. Intended use of the Device .....5
  - C. Additional information .....5
- 2. Safety information .....5
- 3. Technical specifications .....7
- 4. Installation and method for use .....9
  - A. Preliminary precautions .....9
  - B. Installing the fixation system .....10
- 5. Maintenance and after-care .....13
  - 5.A List of authorized cleaning agents and disinfectants .....13
- 6. End of service life .....14
- 7. Product reference .....14
- 8. Applicable standards and directives .....14
- 9. Warranty .....14

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*April 2017***User liability**

This user manual describes how the single-use BiliCocoon™ Fixation System – referred to simply as the “fixation system” (mounting system) – should be installed and used. The fixation system should be used with the BiliCocoon™<sup>1</sup> phototherapy system: the system comprises the fixation system and the phototherapy system, and is also referred to as the “Device” in this document.

The fixation system should only be used appropriately by a trained user and should not be used if it is defective, damaged or if there are any missing parts. In this case, NeoMedLight recommends you dispose of the device and replace it.

The user will be considered solely responsible for the risks incurred by the patient, caregivers, third parties or goods; or if the results differ from those expected due to improper or unsuitable use, or for any improper modifications made to the fixation system.

Any serious incident which occurs in relation to the device should be reported to the manufacturer and the local competent authority of the state in which the user and/or patient is established. You can contact NeoMedLight at the following address: *quality@neomedlight.com*.

According to U.S. Federal Law, the Product may only be sold to authorised professionals and upon their request.

NeoMedLight declares that this product complies with Directive 93/42/EEC.

Year the CE marking was added: 2016.

<sup>1</sup> Also referred to as the BiliCocoon™ System in this document.

# **1. PRODUCT DESCRIPTION AND APPLICATION**

## **A. Description**

The fixation system – or BiliCocoon™ Fixation System – is a fixation system designed for exclusive use with the BiliCocoon™ System, a medical device offering a treatment solution for unconjugated hyperbilirubinemia in new-born babies and which uses phototherapy.

The BiliCocoon™ System consists of a light pad connected to a Light Box emitting blue light in the absorption spectrum of bilirubin, thereby enabling its concentration to be reduced. *For more information, see the BiliCocoon™ System user manual.*

The BiliCocoon™ Fixation System should only be operated by caregivers or a person who has been trained by a qualified trainer.

*For more information, see “3. Technical specifications”.*

## **B. Intended use of the Device**

The BiliCocoon™ Phototherapy System is intended for the treatment of unconjugated hyperbilirubinemia in the population of neonates and infants under 3 months old and weighing less than 10kg. It can be used in the clinical setting or in the home.

## **C. Additional information**

- Only for the treatment of unconjugated hyperbilirubinemia, which is responsible for the majority of cases of neonatal jaundice.
- Only with medical prescription
- Only for the duration indicated in the treatment protocol
- Target patient population: neonates and infants <sup>2</sup>.
- Over the patient's entire body except for genitals and eyes – eye protection must be worn.
- To be used with a single patient
- To be used in hospital or at home: adapted to different environments for use (in an incubator, on a table, in a cradle, in the arms of a parent or caregiver).
- Not to be used on damaged skin

# **2. SAFETY INFORMATION**

---

<sup>2</sup> Under 3 months old and weighing less than 10 kg – definition according to EN 60601-2-50. The term “infant” used in this document is used in compliance with this definition according to EN 60601-2-50.

## Types of safety advice

### **Warning**

Signals danger or risk of injury for the child or the operator

### **Important**

States the necessary actions for obtaining good clinical results

#### **Prerequisites**

The Device must only be used by qualified medical staff who are fully aware of the risks and advantages of this type of treatment and this type of apparatus, or by individuals who have been trained by medical staff.

Warning

**Use at home:** Before using the Device for the first time, it is vital that home users be trained in how to install and use the Device, and be made aware of the risks associated with using it. It is of the utmost importance that they follow the recommendations and instructions provided in the instructions (system and single use) and by healthcare professionals.

The device must not be used if contraindicated, or for medical indications or parts of the body for which it is not recommended or authorised. This could either render the system completely ineffective (no treatment benefits) or could put the patient's health at risk (adverse effects).

Warning

Keep out of reach of children.

Warning

#### **Variation in ambient parameters**

Variation in ambient parameters (room temperature and humidity, exposure to sunlight, surrounding equipment) may have a negative effect on the child's health, especially their temperature or fluid balance.

If the Device is used other than in accordance with instructions, it may malfunction.

Warning

#### **Inflammability**

Keep away from flames or sources of ignition.


**For more information on the risks associated with phototherapy, see the BiliCocoon™ System user manual.**

Warning

#### **Label integrity and readability**

If labels are damaged, you can find information concerning the safety of patients, users, third-parties or property in this document. Contact the supplier if information concerning the identification of the product or of associated products is not visible anymore.

Symbol	Meaning	Symbol	Meaning
	Manufacturer		CE marking declaration according to Directive 93/42/EEC
	See instructions for use		Date of manufacture
	Catalogue reference		Serial number
	Keep dry		Storage and handling temperature range
	Storage and handling humidity range		Fragile, to be handled with care

### **Abbreviations and units of measure**

Abbreviations	Meanings and units of measure
°C	Degrees Celsius (Unit of temperature)
kg	Kilogram (Unit of mass)
cm	Centimetre (Unit of length)
%RH	Humidity (Percentage of humidity in the air)

## **3. TECHNICAL SPECIFICATIONS**

The fixation system – or BiliCocoon™ Fixation System – is a fixation system designed for exclusive use with the BiliCocoon™ System, a medical device offering a treatment solution for unconjugated hyperbilirubinemia neonates and infants and which uses phototherapy.

The BiliCocoon™ System consists of a light pad connected to a Light Box emitting blue light in the absorption spectrum of bilirubin, thereby enabling its concentration to be reduced. *For more information, see the BiliCocoon™ System user manual.*

Parameter	Value	Unit
Dimensions	18 x 15 x 11	cm
Weight	0.6	kg
Colour	Grey	
Min. and max. temperature for use	From 0°C to 40°C	°C
Min. and max. humidity for use	From 30% to 75%	%RH
Min. and max. temperature for storage	From -20°C to 60°C	°C
Min. and max. humidity for storage	From 0% to 90%	%RH

Warning

The disregard of the operating and storage conditions may lead to risks for the patient, the users, third parties or goods.

Warning


Conditions for use and storage are different for the BiliCocoon™ System and the BiliCocoon™ Fixation System. Please ensure that the product is stored and used in accordance with the recommended conditions.


Warning

No not alter the BiliCocoon™ Fixation System. Do not use any accessories, spare parts and equipment not described in the conditions for use. If the fixation system is modified or used improperly, or used with other parts, patients, users, third parties and goods may be at risk.

• **BiliCocoon™ Fixation System**

- 1: Holding claw
- 2: Fixation plate
- 3: Fixation hooks

- 4: Groove
- 5: Fixation screw
- 6: Tightening hinge

## 4. INSTALLATION AND METHOD FOR USE

### A. Preliminary precautions

**Warning**

Before positioning the BiliCocoon on the fixation system, check the Pad's irradiance (see the BiliCocoon™ System user manual).

If the check is not carried out before each use, there is risk of under-exposure, and therefore ineffectiveness of the treatment, or over-exposure.

**Warning**

Ensure that the fixation system is not chipped or damaged, and that there are no sharp edges that could cause harm.

**Warning**

Ensure that the fixation system is clean and unmarked. If necessary, clean it with a suitable cleaning product, such as water or the products specified in this manual.

**Warning**

Ensure that the holding claw is properly attached to the fixation bracket and that there are no missing screws. Check that the mobile hinge can be tightened without any obstruction.

**Warning**

Check that the holding claw is not damaged nor chipped.

**Warning**

Also check that the two fixation hooks are in good condition and are not damaged nor chipped.

**Warning**

The fixation system is designed to be attached to a stand (mounting pole or bed). Check that this stand is solid, stable and sturdy, and that it is between 1 and 5 cm wide. Check the stand (fixation pole or bed) is stable before and after mounting the system. Check the full BiliCocoon™ System is stable before using the system to avoid the fall of the BiliCocoon™ System.

**Warning**

Check that the holding claw is properly screwed on and that it is flush with the stand.

**Warning**

Check that the fixation system is mounted at a suitable height and that it is positioned so it can be easily reached. Check the Light Box™ ventilation holes are not obstructed.

**Warning**

Check that the Light Box™ is correctly positioned on the fixation system.

**Warning**

Check the Light Box™ power lead is not caught on anything.

**Warning**

Failure to follow the conditions of use, maintenance and storage may put patients, users, third parties or goods at risk.

**Warning**

The fixation system – or BiliCocoon™ Fixation System – is designed for exclusive use with the BiliCocoon™ System, a medical device offering a treatment solution for unconjugated hyperbilirubinemia in new-born babies using phototherapy.

**Warning**

While BiliCocoon™ Fixation System is used to support BiliCocoon™ system, the support where the fixation system is fixed should not be moved. Improper use may actually lead to an ineffective (not beneficial) treatment or to risks for the patient's health (unwanted adverse events) and the users.


## **B. Installing the fixation system**

**Warning**


The BiliCocoon™ Light Box can only be installed on the BiliCocoon™ Fixation System once the pads have been removed from the Light Box. The pads should only be connected or disconnected when the BiliCocoon™ Light Box is mounted on the BiliCocoon™ Fixation System. Users must only disconnect or connect the pads while securely holding the Light Box tightly with their hand.

**Presentation of the fixation system**


The fixation system consists of a fixation plate (1) with two fixation hooks, as well as a holding claw (2). Completely unscrew the screw so that the tightening hinge can be fully opened. If the hinge cannot be fully opened, replace the fixation system.

**Hanging the fixation system**

Place the tightening hinge on a stable stand (3), so that the hinge can be tightened on to it. The stand must be stable, solid, vertical and must be between 1 and 5 cm thick. Then, tighten the hinge by fully screwing the tightening screw in as far as it will go. The fixation system must be flush with the stand and must not slide or move.

**Hanging the BiliCocoon™ Light Box**

**Using the BiliCocoon™ Fixation System groove**

Once the fixation system is properly attached to the stand, hang the BiliCocoon™ Light Box (4) by resting the Light Box unit handle on the two fixation hooks. Everything must be stable and well-centred.

The BiliCocoon™ Nest Pad or the BiliCocoon™ Bag Pad (5) can be placed in the groove of the fixation system.

## **5. MAINTENANCE AND AFTER-CARE**

Please clean the BiliCocoon™ Fixation System after each use.

Never immerse in a bath and never spray a solution directly onto the BiliCocoon™ Fixation System's surface.

Do not autoclave or gas sterilise.

If you do not sufficiently clean the device and/or do not follow the cleaning instructions, you risk contaminating or damaging it.

**Warning**

Maintenance should not be carried out on the BiliCocoon™ Fixation System when it is positioned on the fixation system.

Remove the BiliCocoon™ Fixation System from its stand before cleaning it. Use a soft cloth and cleaning and disinfectant agent (see the list of cleaning and disinfectant agents used in paragraph 5.A), then carefully rinse with a soft, wet cloth. Perform a visual inspection of the BiliCocoon™ Fixation System to assure that it is clean. Otherwise repeat the previous instructions until there is no more soil observed. Only use the fixation system once it is completely dry.

Maintenance must not be carried out when the system is connected to the power supply.

### **5.A List of authorized cleaning agents and disinfectants**

The system can be cleaned without being damaged using the following products:

- An EPA-registered cleaning and disinfectant product, such as CaviWipes™, METREX
- A product containing didecyldimethylammonium chloride, such as Surfa'Safe, Laboratoires ANIOS

**Warning**

Do not use products containing:

- Caustic or abrasive cleaning products
- Powerful cleaning solutions (acidic and alkaline)
- Iodine-based solutions
- Phenolic compound-based products, which can lead to increased bilirubin levels in newborn babies

## **6. END OF SERVICE LIFE**

If the fixation system becomes damaged or worn, or no longer works, you should stop using it and replace it; it should be disposed of in accordance with the law in force. The fixation system is made primarily of aluminium.

**Warning**

Do not use the BiliCocoon™ Fixation System if its integrity and/or its functions are compromised by a defect or worn parts.

## **7. PRODUCT REFERENCE**

BCC010	BiliCocoon™ Light Box
BCC101	BiliCocoon™ Nest System
BCC111	BiliCocoon™ Nest Pad
BCC12150	BiliCocoon™ Disposable Nest x50
BCC201	BiliCocoon™ Bag System
BCC211	BiliCocoon™ Bag Pad
BCC22150	BiliCocoon™ Disposable Bag x50
BCC301	BiliCocoon™ Fixation System

## **8. APPLICABLE STANDARDS AND DIRECTIVES**

EN 60601-1	Directive 93/42/EEC	REACH 1907/2006
------------	---------------------	-----------------

## **9. WARRANTY**

The BiliCocoon™ Fixation System comes with a one-year contractual warranty.

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