

# AeroNOx™

Nitric Oxide Delivery

NOx on the go



Compact and portable  
NO and NO<sub>2</sub> monitoring  
Built-in sampling system  
Audible and visual alarms

### Technical Specifications

<b>Standard Features</b>	Light weight, portable design Continuous delivery of inhaled Nitric Oxide Monitors NO, NO <sub>2</sub> , and O <sub>2</sub> Backlit, digital readout Long battery life Push button zeroing Compatible with conventional and Bronchotron® / TXP-2D® high frequency ventilation
<b>Audible and Visual Alarms</b>	NO high / low NO <sub>2</sub> high 1 minute silence feature
<b>Display</b>	Backlit LCD main screen LED indicators
<b>Visual Indicators</b>	NO and NO <sub>2</sub> in ppm O <sub>2</sub> % NO delivery flow rate in lpm Battery indicator Sample line occlusion
<b>Operating Temperature Range</b>	5°C to 40°C
<b>Atmospheric Range</b>	0 m to 3658 m
<b>Dimensions</b>	25.4 cm H x 30.5 cm W x 12.7 cm D
<b>Weight</b>	5.3 kg
<b>Gas Supply</b>	Compatible with standard and portable Nitric Oxide cylinders with CGA 626 fittings
<b>Gases Detected</b>	Nitric Oxide (NO) Nitrogen Dioxide (NO <sub>2</sub> ) Oxygen (O <sub>2</sub> )
<b>Concentration Range</b>	NO 0 to 100 ppm NO <sub>2</sub> 0 to 12 ppm O <sub>2</sub> 5 to 100%
<b>NO Mass Flowmeter</b>	0 to 2 lpm, 10 cc/min resolution
<b>Calibration</b>	Push button zeroing and simple high range calibration
<b>Sensors</b>	Electrochemical
<b>AC Power</b>	Universal AC adaptor, 100-240V AC, 50-60 Hz, 0.75-0.35A
<b>Internal Battery</b>	One, 6V, 4.5-amp-hour sealed, lead acid type with auto recharge 5 hours fully charged

### Worldwide

#### International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754  
T. +1.512.873.0033 F. +1.512.873.9090  
sales@int-bio.com

### Europe

#### International Biomedical AG

Glütschbachstrasse 2, CH-3661 Uetendorf, Switzerland  
T. +41.33.345.66.00 F. +41.33.345.66.08  
welcome@int-bio.ch

### int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.