

Kangaroo Care Phototherapy

The phototherapy system designed for kangaroo care that allows for mother/baby bonding time during treatment.



NICU



L&D

Kangaroo care compatible
Meets and exceeds AAP guidelines
Large surface area minimizes treatment time
Double the surface area for full body coverage



35 $\mu\text{W}/\text{cm}^2/\text{nm}$ intensive irradiance
460 nm peak wavelength
1200 cm^2 (186 in^2) effective surface area

Even light distribution; no hot spots
Soft woven and flexible optical fibers
Quiet operation (< 37 dB(A) at 1 m)



NEST SYSTEM

Ideal for preterm babies
Quick, easy access to newborn
Large surface area with posterior & lateral coverage: 40 x 30 cm
Incubator & warmer compatible



BAG SYSTEM

Ideal for full term babies
Compatible with developmental care
Anterior & posterior pads for full body coverage: 20 x 30 cm (x2)
Residual light minimized

Order Information

Part #	Component	Part #	Component
731-0340	BiliCocoon Nest System	731-0346	BiliCocoon Nest Disposable Pad Covers (50/box)
731-0341	BiliCocoon Bag System	731-0347	BiliCocoon Bag Disposable Pad Covers (50/box)
731-0345	BiliCocoon Mounting System	731-0348	BiliCocoon Spare Parts Kit (10 fuses and filters/box)

Europe

International Biomedical AG

Glütschbachstrasse 100, CH-3661 Uetendorf, Switzerland
T. +41.33.345.66.00 F. +41.33.345.66.08
welcome@int-bio.ch

Worldwide

International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754
T. 512.873.0033 F. 512.873.9090
sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.