



Simplified Neonatal Brain Monitoring





CFM/aEEG specifically designed for NICU clinicians 1, 2 or 8 channels of bedside brain monitoring Simple setup: only two clicks to start monitoring Actively shielded components produce high quality data

BIOMEDICAL

nëo



The Nëo Infant Brain Monitor is an easy-to-use CFM/aEEG solution optimized for every clinician in the NICU. Specifically designed for neonates, the Nëo can provide 1 or 2 channels for routine monitoring or up to 8 channels for more advanced monitoring.

The Nëo's simple touchscreen interface allows clinicians to visualize the aEEG patterns in great detail at the bedside. Powered by a high performance, medical-grade computer and Windows 10 platform, the Nëo provides multiple data export options and advanced data encryption for extra security.

Actively Shielded Components



Neonatal Waveguard[™] Caps The actively-shielded caps provide a non-invasive solution for accurate electrode positioning.



Access To Single Lead Electrodes The color-coded interface cable provides compatibility to hydrogel, goldcup, or needle electrodes.



Research Grade Amplifier The research grade EEG amplifier provides high input impedance for improved signal quality and less artifact.

Diagnostic Support Graphics



The **Inter-Burst-Interval (IBI)** quantifies the duration between bursts of brain activity and can be used to help assess the severity of HIE.



The **Burst-Supression-Ratio (BSR)** indicator represents the fraction of time the brain has spent in suppression and is highly predictive of adverse outcomes following hypoxic-ischemic injury.

Worldwide

International Biomedical 8206 Cross Park Drive, Austin, Texas USA 78754 T. 512.873.0033 F. 512.873.9090 sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.

ISO 13485 CERTIFIED NEO-062020