

TecoTHERM NEO

Infant Cooling System

*Infant cooling from
transport to NICU*



NICU



Ground



Air

Programmed hypothermia protocol
Baby servo controlled
Compact and lightweight
USB data output port

TecoTHERM NEO

Infant Cooling System

Features


Three Operational Modes

- Servo-control complete treatment
- Servo-control constant temperature
- Constant mattress temperature



Roll stand option for convenient and portable total body cooling.

Small Transport Package



*Temperatures based on data from the TOBY Study in the New England Journal of Medicine***

Order Information

Part #	Component	Part #	Component
731-0154	TecoTHERM Neo with Starter Kit*	731-0313*	Reusable Adapter Cable for Skin Temp Probe
731-0255*	Non-Vented Fill-up Set	731-0152*	Reusable Adapter Cable for Esoph/Rectal Probe
731-0501*	Disposable Blankets (3/pack)	731-0240*	Disposable Esophageal/Rectal Probes (20/pack)
731-0147*	Mattress Hose, 2m	731-0312*	Disposable Skin Temp Probes (20/pack)
731-0146	Mattress Hose, 1m	731-0257*	Reusable Cleaning Loop
731-0258	Chlorine Dioxide Tablets (30/pack)	731-0438	Felt Ties (10/pack)

*Parts included in TecoTHERM Neo Starter Kit. Additional parts included in Starter Kit: Power Cable, Chlorine Dioxide Tablets (6), Instruction Manual, and Set-up Guide.

**Moderate Hypothermia to Treat Perinatal Asphyxial Encephalopathy; Azzopardi D, Strohm B, Edwards A, Dyet L, Halliday H, Juszczak E, Kapellou O, Levene M, Marlow N, Porter E, Thoresen M, Whitelaw A, Brocklehurst P; N Engl J Med, 2009; 361:1349-1358.

Europe

International Biomedical AG

Glütschbachstrasse 100, CH-3661 Uetendorf, Switzerland
T. +41.33.345.66.00 F. +41.33.345.66.08
welcome@int-bio.ch

Worldwide

International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754
T. 512.873.0033 F. 512.873.9090
sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.

**ISO 13485
CERTIFIED**

AB008-10-2020