

Positioning Aids

Squishon 2 and 3 and Squishon Mattress



Squishon 2 and 3

Provides an exceptionally soft, pliable surface to support and cradle a baby's head and body. The rectangular shape allows for head movement while maintaining a supportive surface. Made from a soft material and filled with a non-toxic gel.

Squishon 2 and 3

Squishon 2 to support the extremely low birth weight infant. Squishon 3 for older preemies and babies with more resilient skin.

Product #	Legacy #	Description	Quantity per case
1 989803196361	91033-S	Squishon 2 6 x 9 in (15 x 22.5 cm) ¹	12
2 989805606461	91033-3	Squishon 3 6 x 9 in (15 x 22.5 cm) ¹	12
989805606481	91033-A	Disposable covers, Squishon 2, 3 & Gel-E Donuts 7 x 9 in (18 x 23 cm)	50

Squishon Mattress

Pliable, gel-filled mattress that cushions the infant and absorbs vibration.

Product #	Legacy #	Description	Quantity per case
3 989805603301	91017	Squishon Mattress 12 x 22 in (30 x 55 cm) ¹	1
989805606321	91017-A	Disposable covers, Squishon Mattress 12.5 x 28 in (31.8 x 71.1 cm)	25

Squishon Mattress

A gel-filled mattress provides a soft surface designed to help relieve pressure points, protect fragile skin and encourage proper positioning. Caregivers can reposition the baby and shift pressure points without having to handle the infant.

1. Has an 18-month shelf life from the date of manufacture.