


Jaundice Management

WeeSpecs



1

WeeSpecs phototherapy mask

WeeSpecs' two-piece design is made from a soft, hypoallergenic, light-blocking material that can be adjusted for a custom fit. The inside frame of the goggles seals around the skeletal socket, reducing pressure and allowing the eyes to open and close naturally.

WeeSpecs Phototherapy Mask

The WeeSpecs offer protection for baby's eyes with a two-piece, adjustable design.

Product #	Legacy #	Description	Quantity per case
1 989805607631	99071-01	WeeSpecs preemie – fits head circumference < 28 cm (< 11 in)	50
1 989805607641	99071-02	WeeSpecs small – fits head circumference < 28 to 34 cm (< 11 to 13.4 in)	50
1 989805607651	99071-03	WeeSpecs medium – fits head circumference > 34 cm (> 13.4 in)	50

WeeSpecs Supreme Phototherapy Mask

WeeSpecs Supreme feature an easy-to-adjust strap and snug fit.

Product #	Legacy #	Description	Quantity per case
2 989805604201	1022894	WeeSpecs Supreme preemie – fits head circumference < 26 cm (< 10.2 in)	50
2 989805604161	1022481	WeeSpecs Supreme small – fits head circumference < 26 to 34 cm (10.2 to 13.4 in)	50
2 989805604211	1022896	WeeSpecs Supreme medium – fits head circumference > 34 cm (> 13.4 in)	50



2

WeeSpecs Supreme phototherapy mask

WeeSpecs Supreme incorporate an anatomically contoured goggle design and dark fabric that provide superb light-blocking protection without placing pressure on the eyes. WeeSpecs are held in place with a layer of breathable, slip-resistant foam that reduces movement and are designed to fit newborns to approximately three weeks old.