



Infant Resuscitator



NICU



L&D

T-Piece and auxiliary flowmeter
Integrated blender and suction
Meets NRP resuscitation guidelines
All-in-one portable solution



Standard Features



T-Piece or Bag & Mask

Provides an integrated air/oxygen blender to power two delivery options: T-Piece and auxiliary flowmeter.



Complete Resuscitation Package

Resuscitation, blender, and suction all in one. Meets the latest NRP guidelines.

Safety Stops

PIP and flow safety stops prevent distention of the lung with an option to override if higher flow and/or pressure is needed.



Roll stand option for convenient, portable resuscitation.

Order Information

Part #	Component	Part #	Component
738-1099-1	Puffin Infant Resuscitator, DISS	738-1657	Vacuum Filter
738-1099-2	Puffin Infant Resuscitator, NIST	738-1704	Vacuum Canister Ring Mount
738-1099-3	Puffin Infant Resuscitator, Canadian DISS	738-1699	T-Piece Circuit with #1 Size Mask (10/pk)
738-1701	Disposable Vacuum Canister	739-1065	Puffin Roll Stand
738-1702	18" Vacuum Tube	738-1195	A750i Resuscitator Swivel Mount

Europe

International Biomedical AG

Glütschbachstrasse 100, CH-3661 Uetendorf, Switzerland
T. +41.33.345.66.00 F. +41.33.345.66.08
welcome@int-bio.ch

Worldwide

International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754
T. 512.873.0033 F. 512.873.9090
sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.