



Infant Resuscitator

- T-piece resuscitation
- Auxiliary flow meter
- Meets AAP resuscitation guidelines
- Integrated oxygen blender and suction



Everything you need, when you need it the most

Features



Complete Resuscitation Package

Resuscitation, blender, and suction all in one. Meets the latest NRP guidelines.




Safety Stops

PIP and flow safety stops prevent distention of the lung with option to override if higher flow and/or pressure is needed.



T-Piece or Bag & Mask

Equipped with both t-piece and bag and mask outlet. Both outlets are connected to the integrated blender so the proper FiO₂ is always delivered.



Order Information

Part #	Component	Part #	Component
738-1099-1	Puffin Infant Resuscitator, US DISS	738-1657	Vacuum Filter
738-1099-2	Puffin Infant Resuscitator, NIST	738-1704	Vacuum Canister Ring Mount
738-1099-3	Puffin Infant Resuscitator, Canadian DISS	738-1699	T-Piece Circuit with #1 Size Mask (10pk)
738-1701	Disposable Vacuum Canister	739-1065	Puffin Roll Stand
738-1702	45cm Vacuum Tube	738-1195	A750i Resuscitator Swivel Mount