


Delayed Cord Clamping System



T-piece resuscitation while cord is connected*
Warming mattress*
APGAR timer with audible alerts
“V” shaped elevating mattress platform

Features



T-Piece Resuscitation*

T-Piece resuscitation and suction options are available for delayed cord clamping during high risk deliveries. Bedside resuscitation allows for delayed cord clamping and resuscitation to be performed simultaneously.



Elevating Base

Adjustable height mattress platform is designed to suit most types of deliveries. Its ergonomic "V" shape and compact design allows clinician to easily position the LifeStart around the mother.



Storage and Power

A storage basket keeps necessary supplies organized and at the bedside. Two AC outlets provide power for auxiliary support equipment.



Timer

The timer provides distinct, audible alerts at 1, 2, and 3 minutes to indicate how much time has passed with the umbilical cord still intact.



*Optional

Europe

International Biomedical AG

Glütschbachstrasse 100, CH-3661 Uetendorf, Switzerland
T. +41.33.345.66.00 F. +41.33.345.66.08
welcome@int-bio.ch

Worldwide

International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754
T. 512.873.0033 F. 512.873.9090
sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.