



Your passion, our commitment

Caring for mother and baby in labor and delivery,
the NICU, and the well-baby nursery

Table of Contents



Positioning Aids

Developmentally appropriate positioning solutions that support clinicians' efforts in caring for the most fragile infants

1



Calming and Soothing

Innovative, evidence-based solutions designed to comfort even the smallest babies

8



Therapy Support

Simple products that support clinical interventions in the NICU

10



Jaundice Management

Proven jaundice management solutions based on our developmentally friendly approach to treating newborns

15

Positioning Aids

Beanie Bumper



These flexible positioning aids can be formed around the infant to create a weighted boundary and provide support where it is needed. They can provide full containment while still allowing the clinician to visualize the patient.

Beanie Bumper

The Beanie Bumper is universally flexible and comes in two sizes with the choice of reusable or disposable covers.

Product #	Legacy #	Description	Quantity per case
729-0002	N/A	Beanie Bumper reusable positioning aid with a permanent, wipeable urethane cover 40" (101.6cm) long, 4" (10.1cm) diameter	4
729-0003	N/A	Beanie Bumper reusable positioning aid with a permanent, wipeable urethane cover 48" (121.9cm) long, 5" (12.7cm) diameter	4
729-0005	N/A	Disposable covers for Beanie Bumper 40" (101.6cm) long, 4" (10.1cm) diameter	25
729-0007	N/A	Disposable covers for Beanie Bumper 48" (121.9cm) long, 5" (12.7cm) diameter	25
729-0009	N/A	Reusable cloth covers for Beanie Bumper 40" (101.6cm) long, 4" (10.1cm) diameter	25
729-0011	N/A	Reusable cloth covers for Beanie Bumper 48" (121.9cm) long, 5" (12.7cm) diameter	25

Bendy Bumper



These flexible positioning aids mimic the boundaries of the womb to promote containment and flexion, providing a reflex stimulus for extremity extension and recoil.

Bendy Bumper

The Bendy Bumper comes in two sizes and has unique molding capabilities to maintain their form until reshaped by the caregiver.

Product #	Legacy #	Description	Quantity per case
989805604581	1029244	Bendy Bumper complete reusable unit with a permanent, wipeable urethane cover 28 in (71 cm) long, 2.5 in (6.4 cm) diameter	6
989805604571	1029243	Bendy Bumper complete reusable unit with a permanent, wipeable urethane cover 36 in (91.5 cm) long, 4 in (10.2 cm) diameter	6

Bendy Bumper Cover

Two types of covers are available: washable cloth covers for multiple-patient use and disposable covers for single-patient use.

Product #	Legacy #	Description	Quantity per case
989805606761	92065-B2	Reusable cloth covers for Bendy Bumper 28 in (71 cm) long, 2.5 in (6.4 cm) diameter	12
989805606771	92065-B3	Disposable covers for Bendy Bumper 28 in (71 cm) long, 2.5 in (6.4 cm) diameter	50
989805606701	92065-A2	Reusable cloth covers for Bendy Bumper 36 in (91.5 cm) long, 4 in (10.2 cm) diameter	12
989805606711	92065-A3	Disposable covers for Bendy Bumper 36 in (91.5 cm) long, 4 in (10.2 cm) diameter	50

SnuggleUp



The SnuggleUp is a soft, cozy nest that helps provide proper positioning for preterm and ill babies. Its padded footroll and soft, adjustable straps allow movement while gently maintaining appropriate positioning.

Reusable SnuggleUp

SnuggleUps are made from a soft, washable cotton/polyester interlock knit that can be washed and reused.

Product #	Legacy #	Description	Quantity per case
989805606311	91015-X	Extra small for infants < 900 g (< 2 lb)	12
989805606281	91015-A	Small for infants 800 to 1350 g (1.8 to 3 lb)	12
989805606291	91015-B	Medium for infants 1300 to 2250 g (2.8 to 5 lb)	12
989805606301	91015-C	Large for infants > 2250 g (> 5 lb)	12

Disposable SnuggleUp

Disposable SnuggleUps are particularly useful for admissions, procedures, transport, or when laundry facilities are not available.

Product #.	Legacy #	Description	Quantity per case
989805607421	97013-XDIS	Extra small for infants < 900 g (< 2 lb)	12
989805607391	97013-ADIS	Small for infants 800 to 1350 g (1.8 to 3 lb)	12
989805607401	97013-BDIS	Medium for infants 1300 to 2250 g (2.8 to 5 lb)	12
989805607411	97013-CDIS	Large for infants > 2250 g (> 5 lb)	12

Frederick T. Frog



Our versatile, frog-shaped beanbags were designed by a NICU nurse to achieve clinicians' goals for developmentally supportive positioning. Place "Freddy" around the head, neck, or hips, or use it to position an extremity and provide tactile support. Available in two versions: washable/reusable, or wipeable/reusable with launderable or single-use covers.

Frederick T. Frog – Washable

The washable Freddy Frog can be processed through traditional laundry services.

Product #	Legacy #	Description	Quantity per case
989805603341	97040	Frederick T. Frog, washable	6

Frederick T. Frog – Wipeable

The wipeable Freddy Frog does not need to be machine washed, just simply wipe it clean. Choose from a cloth/reusable or a disposable cover.

Product #	Legacy #	Description	Quantity per case
989805635141	1084677	Frederick T. Frog, wipeable with reusable covers	6
989805635151	1084678	Reusable covers	50
989805635131	1084679	Disposable covers	50

Prone Plus



The Prone Plus provides ventral support to infants with a comfortable, hourglass-shaped pad. They support newborns in the proper prone position by using gravity to facilitate natural rounding of the shoulders, hand-to-mouth coordination, and head-and-trunk alignment.

Prone Plus

Provides standard and consistent prone support while helping to maintain a secure, balanced and flexed position.

Product #	Legacy #	Description	Quantity per case
989805604811	1034167	Prone Plus with knit cover, extra small < 800 g (< 1.8 lb)	6
989805604821	1034168	Prone Plus with knit cover, small 800 to 1250 g (1.8 to 2.8 lb)	6
989805604831	1034169	Prone Plus with knit cover, medium 1250 to 2000 g (2.8 to 4.4 lb)	6
989805604891	1034986	Prone Plus with knit cover, large > 2000 g (> 4.4 lb)	6

Prone Plus Covers

Soft, washable knit covers are available for each size of the Prone Plus.

Product #	Legacy #	Description	Quantity per case
989805604841	1034170	Washable knit covers, extra small	6
989805604851	1034171	Washable knit covers, small	6
989805604861	1034172	Washable knit covers, medium	6
989805604901	1034987	Washable knit covers, large	6

Prone Plus Disposable Covers

Soft, disposable covers are available for each size of the Prone Plus.

Product #	Legacy #	Description	Quantity per case
729-0016	N/A	Prone Plus Disposable Cover, extra small	25
729-0017	N/A	Prone Plus Disposable Cover, small	25
729-0018	N/A	Prone Plus Disposable Cover, medium	25
729-0019	N/A	Prone Plus Disposable Cover, large	25

Squishon 2 and 3 and Squishon Mattress



1



2



3

Squishon 2 and 3

Provides an exceptionally soft, pliable surface to support and cradle a baby's head and body. The rectangular shape allows for head movement while maintaining a supportive surface. Made from a soft material and filled with a non-toxic gel.

Squishon 2 and 3

Squishon 2 to support the extremely low birth weight infant. Squishon 3 for older preemies and babies with more resilient skin.

Product #	Legacy #	Description	Quantity per case
1 989803196361	91033-S	Squishon 2 6 x 9 in (15 x 22.5 cm) ¹	12
2 989805606461	91033-3	Squishon 3 6 x 9 in (15 x 22.5 cm) ¹	12
989805606481	91033-A	Disposable covers, Squishon 2, 3 & Gel-E Donuts 7 x 9 in (18 x 23 cm)	50

Squishon Mattress

Pliable, gel-filled mattress that cushions the infant and absorbs vibration.

Product #	Legacy #	Description	Quantity per case
3 989805603301	91017	Squishon Mattress 12 x 22 in (30 x 55 cm) ¹	1
989805606321	91017-A	Disposable covers, Squishon Mattress 12.5 x 28 in (31.8 x 71.1 cm)	25

1. Has an 18-month shelf life from the date of manufacture.

Squishon Mattress

A gel-filled mattress provides a soft surface designed to help relieve pressure points, protect fragile skin and encourage proper positioning. Caregivers can reposition the baby and shift pressure points without having to handle the infant.

Gel-E Donut



An exceptionally soft skin combined with a unique gel formulation makes Gel-E Donuts soft, yet supportive. They are designed to help relieve pressure caused by prolonged immobility and other conditions where frequent handling may be contraindicated. Reengineered in 2015 to change to a non-porous skin.

Gel-E Donuts

Comfortable & pliable to support and cradle the baby's head.

Product #	Legacy #	Description	Quantity per case
989803196371	92025-XS	Extra small 7 in diameter (18 cm), fill weight: 325 g (11.5 oz), purple ¹	12
989803196381	92025-S	Small 7 in diameter (18 cm), fill weight: 400 g (14 oz), green ¹	12
989803196391	92025-M	Medium 7 in diameter (18 cm), fill weight: 450 g (16 oz), peach ¹	12
989805606481	91033-A	Disposable covers, Squishon 2, 3, & Gel-E Donuts	50

1. Has an 18-month shelf life from the date of manufacture.

Calming & Soothing

Sweet-Ease

This 24% sucrose and purified water solution is a proven way to calm and soothe infants in the NICU, PICU, or well-baby nursery. It's available in a 15 mL cup, 2 mL vial, or 1 mL vial.



1

Sweet-Ease 15 mL cup

This 15 mL cup is aseptically packaged with a peel-off lid that is suitable for dipping a pacifier or for administration via a dropper or syringe.

Sweet-Ease 15 mL cup

Aseptically packed cups give clinicians flexibility for administration of sucrose solution.

Product #	Legacy #	Description	Quantity per case
1 989805603401	99044	Sweet-Ease natural (15 mL cup) ¹	200

Sweet-Ease 1 & 2 mL vials

Conveniently packaged vials are great for quick delivery of sucrose solution to the infant.

Product #	Legacy #	Description	Quantity per case
2 989805637451	1093977	Sweet-Ease natural (2 mL vial) ²	250
2 00850018561235	N/A	Sweet-Ease natural (1 mL vial) ²	250

1. Has a one-year shelf life from the date of manufacture.
 2. Has a two-year shelf life from the date of manufacture.



2

Sweet-Ease 1 & 2 mL vials

Vials simplify clinical workflow by allowing the administration of sucrose directly on to the tip of an infant's tongue or buccal surface.

Snoedel



The Snoedel, originating from a Dutch family tradition, offers parents a means to comfort and bond with their child even when the baby is in the hospital. The pure sheep wool incorporated into the high-quality flannel doll absorbs the parent's scent and slowly releases it, supporting the sensitive olfactory system of the infant.

Snoedel

The Snoedel is a flannel bonding doll and positioning aid.

Product #	Legacy #	Description	Quantity per case
989805607221	96029-B	Snoedel (unisex print)	12

Therapy Support

Infant Heel Warmer



This easy-to-use heel warmer provides steady warming without skin adhesives or wet wash towels. The temperature of this heel warmer is dependable during each use providing consistency for the heel stick procedure.

Infant Heel Warmer

Provides consistent, uniform heat distribution for increased blood flow.

Product #	Legacy #	Description	Quantity per case
989805603201	1223	Infant Heel Warmer	100 (25 per box; 4 boxes per case)

Heel Snuggler



The Heel Snuggler is the only infant heel warmer designed to conform truly to the infant's heel while staying firmly in place without skin adhesives. The unique design is shaped like a butterfly with two soft perforated straps for a comfortable and secure fit.

Heel Snuggler

Unique shape conforms to the infant's heel, keeping warmth where it is needed.

Product #	Legacy #	Description	Quantity per case
989805603411	99047	Heel Snuggler	100 (25 per box; 4 boxes per case)

Infatherm Disposable Warming Mattress



Warming gel mattress designed to provide gentle thermoregulation support for infants big and small. The portable, non-electric heat source provides instant short term warmth for infants during transport or in the delivery room.

Infatherm Disposable Warming Mattress

The durable construction is designed to prevent leaks, while the soft, non-slip fabric protects the infant's delicate skin.

Product #	Legacy #	Description	Quantity per case
989805616831	1015	Infatherm Disposable Warming Mattress	6

Wee Sight Transilluminator



The Wee Sight transilluminator helps to locate the tiny veins in a neonate's small limbs for improved insertion of an IV.

Wee Sight Transilluminator

The Wee Sight uses two AA alkaline batteries (included) and outputs a light wave length of approximately 629 nm.

Product #	Legacy #	Description	Quantity per case
989805603761	1017920	Wee Sight Transilluminator	1

Temperature Trend Indicators



Temperature Trend Indicators provide non-invasive, rapid temperature data for consistent monitoring of a patient's temperature.

Temperature Trend Indicators

Temperature Trend Indicators are accurate to $\pm 1^\circ\text{F} / 0.5^\circ\text{C}$, and are ideal for use in a number of care settings.

Product #	Legacy #	Description	Quantity per case
989805623001	LC105D	LC Accu-Bar Temperature Trend Indicator and dispenser box (one degree thermometer) 94°F-105°F (34.5°-40.5°C)	100

Sali-Wipes



Sali-Wipes sterile saline (isotonic) wipes come in convenient foil packs. Nursing organizations recommend use of sterile isotonic saline for cleansing, as alcohol wipes can be damaging to fragile infant skin and may be absorbed into the bloodstream.

Sali-Wipes

Sali-Wipes are used to clean skin after disinfection & prior to invasive procedures.

Product #	Legacy #	Description	Quantity per case
989805607531	98098-013	Sali-Wipes ¹ 4 x 6 in (10 x 15.2 cm)	576

1. Has an two-year shelf life from the date of manufacture.

BBG Nasal Aspirator



The BBG Nasal Aspirator allows efficient, non-invasive suction of the nasal cavity. It is a safe, easy-to-use bulb-style aspirator that features an over-insertion guard. Two flexible sizes enable caregivers to provide safe suction in a variety of care environments.


BBG Nasal Aspirator

The BBG Nasal Aspirator facilitates gentle suction to remove excess mucus or secretions from a baby's nasal passage.

Product #	Legacy #	Description	Quantity per case
989805606951	93032-02	BBG Nasal Aspirator 3.25 in (8.25 cm)	50
989805606961	93032-S	BBG Nasal Aspirator 4.5 in (11.43 cm)	50

Jaundice Management

WeeSpecs



WeeSpecs Phototherapy Mask

WeeSpecs' two-piece design is made from a soft, hypoallergenic, light-blocking material that can be adjusted for a custom fit. The inside frame of the goggles seals around the skeletal socket, reducing pressure and allowing the eyes to open and close naturally.

WeeSpecs Phototherapy Mask

The WeeSpecs offer protection for baby's eyes with a two-piece, adjustable design.

Product #	Legacy #	Description	Quantity per case
1 989805607631	99071-01	WeeSpecs preemie – fits head circumference < 28 cm (< 11 in)	50
1 989805607641	99071-02	WeeSpecs small – fits head circumference < 28 to 34 cm (< 11 to 13.4 in)	50
1 989805607651	99071-03	WeeSpecs medium – fits head circumference > 34 cm (> 13.4 in)	50

WeeSpecs Supreme Phototherapy Mask

WeeSpecs Supreme feature an easy-to-adjust strap and snug fit.

Product #	Legacy #	Description	Quantity per case
2 989805604201	1022894	WeeSpecs Supreme preemie – fits head circumference < 26 cm (< 10.2 in)	50
2 989805604161	1022481	WeeSpecs Supreme small – fits head circumference < 26 to 34 cm (10.2 to 13.4 in)	50
2 989805604211	1022896	WeeSpecs Supreme medium – fits head circumference > 34 cm (> 13.4 in)	50



WeeSpecs Supreme Phototherapy Mask

WeeSpecs Supreme incorporate an anatomically contoured goggle design and dark fabric that provide superb light-blocking protection without placing pressure on the eyes. WeeSpecs are held in place with a layer of breathable, slip-resistant foam that reduces movement and are designed to fit newborns to approximately three weeks old.

Bili Pod Phototherapy Swaddle



The Bili Pod is a swaddle used during phototherapy treatment with nearly 95% of light transmission through the mesh fabric.

Bili Pod

The Bili Pod meets the needs of babies undergoing phototherapy, withdrawals (NAS), or cooling therapy.

Product #	Legacy #	Description	Quantity per case
731-0360-1	N/A	Bili Pod, Micro Preemie, 1.5-3 lbs	10
731-0360-2	N/A	Bili Pod, Preemie, 3-5 lbs	10
731-0360-3	N/A	Bili Pod, Newborn, 5-13 lbs	10

MRI Compatible Bili Pod

A Bili Pod designed specifically to keep the infant swaddled during a MRI procedure.

Product #	Legacy #	Description	Quantity per case
731-0567-1	N/A	MRI Compatible Bili Pod, Micro Preemie, 1.5-3 lbs	10
731-0567-2	N/A	MRI Compatible Bili Pod, Preemie, 3-5 lbs	10
731-0567-3	N/A	MRI Compatible Bili Pod, Newborn, 5-13 lbs	10



 international
BIOMEDICAL®

For more information visit www.int-bio.com or email sales@int-bio.com.

Worldwide


International Biomedical
8206 Cross Park Drive
Austin, TX 78754 USA

512.873.0033
sales@int-bio.com

Europe

Glütschbachstrasse 100
CH-3661
Uetendorf, Switzerland

+41.33.345.66.00
welcome@int-bio.ch



International Biomedical
8206 Cross Park Drive
Austin, TX 78754 USA

© 2023 International Biomedical. All rights are reserved.
International Biomedical reserves the right to make changes in specifications and/
or to discontinue any product at any time without notice or obligation and will not be
liable for any consequences resulting from the use of this publication. Trademarks
are the property of International Biomedical or their respective owners.

CHMV-3-23

Made in the USA



www.int-bio.com