

# AirBORNE™

750i Transport Incubator

## Neonatal Transport Incubator



Compact and lightweight  
Optional integrated pulse oximetry and O<sub>2</sub> analyzer  
Rugged and reliable

### Technical Specifications

<b>Standard Features</b>	Large double-wall chamber with 360° visibility Head and front access doors 2 hand ports on front door 2 hand ports on back of infant chamber 5 access ports for IV/respiratory tubing Slide out intubation tray Quartz halogen observation light Extended-life battery	Skin temperature probe Large LED display Humidification reservoir IV pole 2 pairs of infant positioning straps Power cord wrap Cylinder mounts Pressure diffusing mattress
<b>Audible and Visual Alarms</b>	System failure Temp sensor failure Air flow failure AC failure	Power failure High temp Low battery
<b>Visual Indicators</b>	Set air temp Set temp >37°C Air temp Baby temp Battery life	Battery power Battery charging AC power External DC power
<b>Dimensions</b>	50 cm H x 55 cm W x 127 cm L (w/o cart, IV pole, or exam light)	
<b>Weight</b>	43 kg (w/o cart, IV pole, or exam light)	
<b>AC Power</b>	230V AC, 50-400 Hz	
<b>External DC Power</b>	12V DC input, 10A maximum	
<b>Internal Battery</b>	One 12V 26-amp hour sealed lead acid gel type with automatic recharge	
<b>Digital Air Temperature Control</b>	17.0°C to 38.9°C range, 0.1°C increments	Display resolution 0.1°C
<b>Sound Level</b>	<60 dB	

### Additional Options Available

<b>Integrated Pulse Oximetry</b>	Masimo® or Nellcor™, and O <sub>2</sub> analyzer with audible and visual alarms	
<b>Conventional Ventilator</b>	Neonatal transport ventilator with remote vent outputs	
<b>Electric Suction</b>	Portable, battery-powered suction	
<b>Patient Monitor</b>	Install hospital-supplied monitor or one provided by International Biomedical	
<b>Thermal / Privacy Cover</b>	Provides patient privacy	Protects from external influences
<b>Resuscitation</b>	Disposable resuscitation bag and mask	T-piece resuscitation
<b>Lighting</b>	LED light bar (430 Lux average)	
<b>Shelves / Mounts</b>	Mounting posts 2.5 cm diameter x 12.7 cm, or 2.5 cm x 19.0 cm	Overhead shelf and / or side platforms
<b>Mobile Cart</b>	2 storage shelves Enclosed bottom shelf Large 5 inch locking casters 4 shock absorbers	

### Worldwide

#### International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754  
T. 512.873.0033 F. 512.873.9090  
sales@int-bio.com

### Europe

#### International Biomedical AG

Glütschbachstrasse 100, CH-3661 Uetendorf, Switzerland  
T. +41.33.345.66.00 F. +41.33.345.66.08  
welcome@int-bio.ch

### int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.