

Infant Resuscitator

- T-piece resuscitation
- Auxiliary flow meter
- Meets AAP resuscitation guidelines
- Integrated oxygen blender and suction



Everything you need, when you need it the most

Features



Puffin®
Infant Resuscitator



Complete Resuscitation Package

Resuscitation, blender, and suction all in one. Meets the latest NRP guidelines.




Safety Stops

PIP and flow safety stops prevent distention of the lung with option to override if higher flow and/or pressure is needed.



T-Piece or Bag & Mask

Equipped with both t-piece and bag and mask outlet. Both outlets are connected to the integrated blender so the proper FiO_2 is always delivered.



Order Information

| Part # | Component | Part # | Component |
|------------|---|----------|--|
| 738-1099-1 | Puffin Infant Resuscitator, US DISS | 738-1657 | Vacuum Filter |
| 738-1099-2 | Puffin Infant Resuscitator, NIST | 738-1704 | Vacuum Canister Ring Mount |
| 738-1099-3 | Puffin Infant Resuscitator, Canadian DISS | 738-1699 | T-Piece Circuit with #1 Size Mask (10pk) |
| 738-1701 | Disposable Vacuum Canister | 739-1065 | Puffin Roll Stand |
| 738-1702 | 45cm Vacuum Tube | 738-1195 | A750i Resuscitator Swivel Mount |

Europe

International Biomedical AG

Glütschbachstrasse 100, CH-3661 Uetendorf, Switzerland
T. +41.33.345.66.00 F. +41.33.345.66.08
welcome@int-bio.ch

Worldwide

International Biomedical

8206 Cross Park Drive, Austin, Texas USA 78754
T. 512.873.0033 F. 512.873.9090
sales@int-bio.com

int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.

ISO 13485
CERTIFIED

LB002-10/2019