

## Kangaroo Care Phototherapy

*The phototherapy system designed for kangaroo care that allows for mother/baby bonding time during treatment.*



NICU



L&D

**Kangaroo care compatible**  
**Meets and exceeds AAP guidelines**  
**Large surface area minimizes treatment time**  
**Double the surface area for full body coverage**



35  $\mu\text{W}/\text{cm}^2/\text{nm}$  intensive irradiance  
460 nm peak wavelength  
1200  $\text{cm}^2$  (186  $\text{in}^2$ ) effective surface area

Even light distribution; no hot spots  
Soft woven and flexible optical fibers  
Quiet operation (< 37 dB(A) at 1 m)



### NEST SYSTEM

Ideal for preterm babies  
Quick, easy access to newborn  
Large surface area with posterior & lateral coverage: 40 x 30 cm  
Incubator & warmer compatible



### BAG SYSTEM

Ideal for full term babies  
Compatible with developmental care  
Anterior & posterior pads for full body coverage: 20 x 30 cm (x2)  
Residual light minimized

## Order Information

Part #	Component	Part #	Component
731-0340	BiliCocoon Nest System	731-0346	BiliCocoon Nest Disposable Pad Covers (50/box)
731-0341	BiliCocoon Bag System	731-0347	BiliCocoon Bag Disposable Pad Covers (50/box)
731-0345	BiliCocoon Mounting System	731-0348	BiliCocoon Spare Parts Kit (10 fuses and filters/box)

### Distributed by

International Biomedical  
8206 Cross Park Drive, Austin, Texas USA 78754  
T. 512.873.0033 F. 512.873.9090  
sales@int-bio.com

### int-bio.com

International Biomedical reserves the right to make changes in specifications and features shown herein or discontinue the product described at any time without notice or obligation. Contact your representative for the most current information.



ISO 13485  
CERTIFIED

AB008-10/2019